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LIABILITY

The user of Like-A implants has the duty to determine whether or not any product is suitable for the particular patient and circumstances. DEEP Company disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any

errors in professional judgment or practice in the use of DEEP Company products. The clinician is also obliged to study the latest developments in regard to the DEEP Company products and their applications. In cases of doubt, the clinician has to contact DEEP Company. Since the processing and surgical application of this

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DEEP Company
 ZA du bois Saint-Pierre
 38280 JANNEYRIAS - FRANCE

Tel : +33 (0) 4 28 01 26 03
www.deep-company.com
contact@deep-company.com

IMPORTANT : Please read these instructions carefully before use

product are under the control of the clinician, they are his/her responsibility. DEEP Company assumes no liability whatsoever for damages arising thereof.

DESCRIPTION

Endosseous threaded dental implant made from biocompatible pure grade 4 titanium.

INDICATIONS

Dental implants are used for replacing teeth both in mandible and jawbone. They must be placed in the jaw bone to be integrated (osteointegration). Restorations range from a single replacement to the use of elements for overdenture applications as well as the full arch of bridgework.

CONTRAINDICATIONS

It is contraindicated placing dental implants in patients:

- medically unfit for an oral surgical procedure
- with inadequate bone volume, unless an augmentation procedure can be considered
- for whom it is impossible to place the adequate size or number of implants, or in the desirable position; treatment plan is then incompatible with the physiological functional constraints.

INSTRUCTION FOR CLINICIAN

This product must be manipulated/placed by trained and qualified persons having taken knowledge of this note.

The surgery must be performed by a clinician mastering the specific surgical techniques for dental placement; we strongly recommend that clinicians always go through special training before undertaking a new treatment method.

The clinician is responsible for any complications that can occurred from a use that do not comply with the indications and contraindications, from a defective surgical technique or a lack of asepsis. These complications shall in no way be attributed to DEEP Company.

Further information on the indication and use of different systems are available at DEEP Company or on www.deep-company.com or www.generictimplants.com.

GENERAL PRECAUTIONS

One hundred percent implant success can never be guaranteed.

With respect to children, the implant placement is not recommended until the end of the jaw bone growth has not been duly noted. Pre-operative deficits in gingival and bone level, may yield a compromised aesthetic result or unfavourable implant angulation.

PRE-OPERATIVE PRECAUTIONS

A careful clinical and radiological examination of the patient has to be performed prior to surgery to determine the physical and psychological status of the patient. Special attention has to be given to patients who have local or systemic factors that could interfere with the healing process of either bone or soft tissue or the osteointegration process (cigarette smoking, uncontrolled diabetes, oro-facial radiotherapy, steroid therapy, infections in the neighboring bone). In general, implant placement and prosthetic design must accommodate individual patient conditions. In case of bruxism or unfavorable jaw relationships, reappraisal of the treatment option may be considered.

PER-OPERATIVE PRECAUTIONS

Particular caution should be used when placing small diameter implant (Ø3.5 NP) in the posterior region due to the risk of prosthetic overload. Lack of adequate quantity and/or quality of remaining bone, infection and generalized diseases may be potential causes for failure of osteointegration both immediately or after the surgery. Besides the mandatory precautions for any surgery (such as asepsis), during drilling in the jaw bone, the clinician must avoid damage the nerves and vessels by referring to anatomical knowledge and preoperative radiographs. Failure to recognize actual lengths of drills relative to radiographic measurement or drilling beyond the depth intended can result in permanent injury to nerves or other vital structures. It can result in a permanent numbness to the lower lip and chin or lead to hemorrhage such as in the floor of the mouth.

Implants may be tilted up to 45° relative to the occlusal plane. When used with angulations between 30° and 45°, the following applies: The tilted implants must be splinted; a minimum of 4 implants must be used when supporting a fixed prosthesis in a fully edentulous arch.

The practitioner is advised to provide regular patient follow up and to inform about good oral hygiene.

PATIENT INFORMATION

The clinician is responsible for:

- informing the patient of the potential risks and adverse effects associated with the fitting of such an implant system,
- informing the patient that the safety and effectiveness of the devices depends on his behaviour and his oral hygiene,
- ensuring that the patient reports him all suspicious movements between implant and abutment, and will have to undergo periodic inspections for detecting signs that precede abnormal functional symptoms.

POTENTIAL ADVERSE EVENTS

Treatment by means of implants may lead to loss of bone, biologic or mechanical failures including fatigue fracture of implants.

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INTERACTION WITH MAGNETIC RESONANCE IMAGING AND OTHER ELECTROMAGNETIC SYSTEM

The raised temperature after interaction between MRI and implants does not pose any known hazard to the patient.

But the quality of the MR image may be compromised, if the area of interest is in the same area of the device: an artefact may appear on the picture.

STERILITY

All implants are delivered sterile and must be used prior to the labelled expiration date. They are for a single use only. Do not use implants if the packaging has been damaged or previously opened. Never reprocess an implant. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

PRODUCT ELIMINATION

Sorting between infectious waste and waste comparable to household waste lies with the professional who produced the waste.

- Typically common waste products such as packaging are similar to household waste and follow the same waste management system.
- Items and materials soiled by blood and / or body fluids have to be treated as waste from riskily treatment activities.

RETURN

For any returned product to the aftersales service, please clean and sterilize the product following described actions before return them to us under a double blister.

SURGICAL PROCEDURES

1. During drilling procedures, bone quality should be considered (see figure A). Drilling must proceed at high speed (maximum 2000 rpm for step/twist drills) under constant and profuse irrigation by sterile saline at room temperature.
2. Prepare implant site. When using a flapless surgery, add-on soft tissue height to drill depth.
3. Measure the final depth of implant site for applicable implant length using a depth probe. Use the drills graduations to perform a drilling with the proper depth.
Caution: Take into account the additional drill length when drilling near vital structures.
4. Open the implant package and pick up the implant from inner casing with the implant driver. The implants are ideally installed with low speed, max 25 rpm, using a surgical machine or a manual torque wrench.
5. Place and tighten the implant without exceeding an insertion torque of 70 Ncm.
Caution: Never exceed the insertion torque of 70 Ncm. Overtightening an implant may lead to damage of the implant, fracture or necrosis of the bone site.
6. If the implant gets stuck during implant installation or 70 Ncm is achieved before fully seated, rotate the implant counter clockwise using drilling machine or manual torque wrench, in reverse mode, and remove the implant from the site. Replace the implant back into titanium casing. Use a wider drill, countersink or screwtap to enlarge the site. If screwtap is used, place it into the implantation site, select the low speed mode (25 rpm) and drill to the appropriate length. Switch to reverse mode and pull the screwtap out. Continue the implant installation until desired position is achieved.
7. For the immediate function, the implant should be able to withstand a final torque of 25-70 Ncm.
8. Depending on surgical protocol, place a cover screw or an abutment and suture.

For any additional information on surgical procedures please consult the "Like-A product documentations" treatment guidelines available at www.generictimplants.com.

A. Recommended drill sequences based on bone quality.

IMPLANT DIAMETER	SOFT BONE TYPE IV	MEDIUM BONE TYPE II-III	DENSE BONE TYPE I
Ø 3.5	2.0 (2.4/2.8) (Tap 3.5)	2.0 2.4/2.8 (2.8/3.2) Tap 3.5	2.0 2.4/2.8 (2.8/3.2) Dense bone Tap 3.5
	Ø 4.3	2.0 2.4/2.8 (2.8/3.2) (Tap 4.3)	2.0 2.4/2.8 (3.2/3.6) (3.8/4.2) Dense bone Tap 4.3
Ø 5.0	2.0 2.4/2.8 3.2/3.6 (Tap 5.0)	2.0 2.4/2.8 (3.2/3.6) Tap 5.0	2.0 2.4/2.8 (3.2/3.6) (3.8/4.2) (4.2/4.6) Dense bone Tap 5.0

Note: all data is stated in mm.

Drills within brackets (-) denote widening of the cortex only, not drilling the full drilling depth.

SYMBOLS	Device manufacturer	Batch number of the Device	Device reference	Have to not be used if the package has been damaged or opened
	Date after which the medical device is not to be used	Do not resterilized	Single use only	CE Mark with the number of the Notify Body 0123
	Consult the instructions for use for important cautionary information	Need protection from light sources	Have to be protected from moisture	sterilized using irradiation

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DEEP Company thanks you for your support and remains at your disposal for any additional information you may require

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